



mini **UNIKO F**

Art. MUN.F

Mariotti introduce la versione economica MiniUniko. Mantenendo le solite elevate performance, semplicità di utilizzo e praticità. Grazie ad un innovativo software MiniUniko.F è programmabile in pochi secondi e dispone di grande flessibilità per l'utilizzo coi vari sistemi implantari.

#### CARATTERISTICHE TECNICHE:

Tensione di alimentazione: 230/125 Vac - 50/60Hz

Potenza assorbita: 109 VA

Torque max: limitato a 60Ncm

Isolamento: classe I, tipo BF

Velocità motore: 200-30000 rpm

Portata max pompa peristaltica d'irrigazione: 90 ml/min

Dimensioni unità di controllo: mm 245x245x100

Pedale: of/off motore o motore+pompa

Funzioni principali: regolazione velocità, torque, portata pompa, rapporto di riduzione manopola, reverse, 4 programmi, aggiornamento software, motore/cavo/connettore autoclavabili

#### ACCESSORI:

Art. C20 Contrangolo per implantologia della migliore qualità. Micro testina, riduzione 20:1, irrigazione esterna e possibilità di irrigazione interna.

Mariotti introduces the MiniUniko economical version. Always with high performance, it combines ease of use and practicality. Miniuniko.F can be set in few seconds thanks to an innovative software and the great flexibility allows to use and adapt the various implantology systems.

#### SPECIFICATIONS:

Power supply: 230/125 Vac - 50/60Hz

Absorbed power: 109 VA

Maximum Torque: limited a 60Ncm

Insulation: Class 1, type BF

Motor speed: 200-30000 rpm

Peristaltic pump max irrigation: 90 ml/min

Control box dimension: mm 245x245x100

Pedal: motor or motor+pump on/off

Main functions: control of speed, torque, peristaltic pump, handpiece reduction, reverse, 4 programs, software upgraded, autoclavable motor/cable/connector.

#### ACCESSORIES:

Art. C20 Implantology contra-angle for finest quality. Micro-head, 20:1 reduction, external sprayer with possible internal spraying.



OPTIONAL

Art. PED

Pedaliera multifunzionale: Attiva il motore, inversione di rotazione, irrigazione on/off, richiamo dei programmi

Multi-function foot-control: Used to start motor, reverse function, on/off spraying and retrieve programs



MARIOTTI & C. srl  
Micromotors & Milling-Machines  
Via Seganti, 73 - 47100 Forlì - Italy  
Tel +39 0543 474105 - Fax +39 0543 781811  
info@mariotti-italy.com  
www.mariotti-italy.com

Member of  
**UNIDI**  
THE ITALIAN DENTAL  
INDUSTRIES ASSOCIATION



**CE**  
0476